



# LABORATORY DRYING OVEN 500 LITER

**AGGREGATE** 

www.testmak.com





# **Laboratory Drying Oven 500 Liter Capacity**



# **PRODUCT MODEL**

| AG1300/N500     | Laboratory Drying Oven 500 Liters Capacity, Natural Convection - 220 V 50/60 Hz  Laboratory Drying Oven 500 Liters Capacity, Natural Convection - 110 V 60 Hz |  |
|-----------------|---|--|
| AG1300/N500-110 |   |  |
| AG1300/F500     | Laboratory Drying Oven 500 Liters Capacity, Forced Convection - 220 V 50/60 Hz  |  |
| AG1300/F500-110 | Laboratory Drying Oven 500 Liters Capacity, Forced Convection - 110 V 60 Hz   |  |

## **PRODUCT STANDARDS**

| Standards | EN 932-5, 1097-5   ASTM C127, C136, D558, D559, D560, D698, D1557, D1559   BS 1377:1, 1924:11 |            |
|-----------|---|------------|
|           | Stanuarus   | UNE 103300 |

# **INFORMATION**

| Manufacturer      | TESTMAK INS.LAB.MAK.SAN.VE TİC. PAZ. ITH. IHR. LTD. STI |
|-------------------|---|
| Country of Origin | TURKEY  |
| Product name      | Laboratory Drying Oven 500 Liters Capacity              |



Telephone : +90 312 395 36 42 Website : www.testmak.com Koca Sinar Fax : +90 312 395 36 01 E-Mail : info@testmak.com 40 06370 0





# **Laboratory Drying Oven 500 Liter Capacity**

### **DESCRIPTION**

The AG1300/500 series Laboratory Ovens have been designed for drying asphalt, soil, rock, concrete, aggregate or similar materials. The oven capacity is 500 liters. From ambient to 250°C temperature range with a precision of ±2 °C. The interior is manufactured from stainless steel and the exterior is robustly constructed from sheet steel finished in powder coated paint. Laboratory Ovens are fitted with direct reading digital control unit and equipped with an analogue over temperature protection thermostat. 500 liters laboratory ovens are supplied complete with 4 shelves.

#### **Natural Convection Model**

Natural convection ovens rely on temperature differences within the oven to transfer heat to samples. They are used for drying applications requiring a gentle airflow in baking, conditioning, curing, pre-heating, and aging. Natural convection ovens are ideal for sterilization, drying, and thermal storage tasks that do not require high drying rates or special time parameters. Look for safety features that prevent overheating, excellent temperature uniformity, fast heat up and recovery times. Additional features that are nice to have are viewing windows, adjustable shelves, door alarms, lockable doors, stoving and curing options and moisture extraction options.

#### **Forced Convection Model**

Forced air convection ovens are built with a fan inside the wall of the oven which forces the hot air in the oven to circulate throughout it. The forced air circulation system provides both exceptional temperature uniformity and rapid heat recovery. They usually operate at temperatures in a range from ambient to a maximum of between 200 and 300 degrees Celsius. Forced-air laboratory ovens are often used in procedures involving cultures and samples when accurately controlled uniform heating is required. Choose from tabletop, undercounter or upright depending on lab space available. Look for convenient features such as adjustable air intake and exhaust vents, sliding window panel covers and adjustable shelves and safety features such as independent overtemperature protection and isolated heating elements.

## **TECHNICAL SPECIFICATIONS**

| Internal Dimensions | 900x750x750 mm     |
|---------------------|--------------------|
| External Dimensions | 1250x760x1100 mm   |
| Shelves             | 4 pcs              |
| Weight (approx)     | 115 kg             |
| Power               | 220-240 V 50/60 Hz |

Telephone: +90 312 395 36 42 Fax: +90 312 395 36 01 Website: www.testmak.com E-Mail: info@testmak.com





Telephone: +90 312 395 36 42 Fax: +90 312 395 36 01 Website: www.testmak.com E-Mail: info@testmak.com Koca Sinan Industrial Sites 1183 Street No: 40 06370 OSTIM / ANKARA / TURKEY